

PRODUCT LINE EXTENSION

NEW STRAIN NOW AVAILABLE

EZ-Accu Shot™

Background

Our #1 Growth Promotion Testing product, EZ-Accu Shot™ continues to grow as we expand the strains offered based on customer requests/feedback. *Penicillium chrysogenum*, a fungus with the potential to contaminate pharmaceutical laboratories, is joining the Microbiologics product family in the EZ-Accu Shot™ format.

Product Details

New Strains:

- 0178A *Penicillium chrysogenum* derived from ATCC® 10106™

Overview:

- Accredited reference material under ISO 17034 standard
- Delivers 10-100 CFU per 0.1 ml as required by USP/Ph. Eur./JP
- Strains are less than or equal to three passages from reference culture – meets USP/Ph. Eur./JP requirements
- Instant dissolve pellet reduces preparation time
- Pending USP chapter requires this strain for Respiration based rapid microbial methods

Kit Components:

- 5 vials of a single, quantitated microorganism (1 lyophilized pellet per vial)
- 5 vials of Hydrating Fluid (1.2 ml in each vial)
- Up to 50 tests per kit
- Illustrated Instructions

PRODUCT LINE EXTENSION

Marketing and Promotions

- Find a variety of Accu Shot and Growth Promotion Testing Flyers, Technical Publications and Videos here:
 - [Promotional Literature](#)
 - [Technical Publications](#)
 - [Growth Promotion Test Guide for Media Used in Tests for Burkholderia cepacia complex](#)
 - [Growth Promotion Test Guide for Specified Organisms](#)
 - [Growth Promotion Test Guide for Sterility Tests](#)
 - [Growth Promotion Test Guide for Media Used in Microbial Enumeration Tests](#)
 - [Video](#)